

# NextDent 5100

## Professional 3D Printer

### Release Notes

Firmware version: 1.1.2

Release Date: April 2021



## 1 INTRODUCTION

We are proud to announce the new printer firmware version 1.1.2!

To download this update, open 3D Sprint and ensure that you have connected to your NextDent 5100 printer. 3D Sprint will check for new firmware and inform you that version 1.1.2 is available. Follow the on-screen instructions to send the update to the printer. You will need to perform these steps for each printer you have.



**NOTE: If you have not previously upgraded to firmware 1.1.1, please follow the [instructions here to upgrade](#), as 3D Sprint did not upgrade the printer directly prior to the 1.1.1 release. Release notes are [linked here](#).**

## 2 NEW FEATURES

### BUTTON IMPROVEMENTS

Improvements have been made to the buttons seen on screen after a print. This improvement has been applied to all languages supported in the firmware.



**3D Systems, Inc.** 333 Three D Systems Circle Rock Hill, SC 29730  
[www.3dsystems.com](http://www.3dsystems.com)

© 2021 by 3D Systems, Inc. All rights reserved. Specifications subject to change without notice. 3D Systems, the 3D Systems logo, ProX and Duraform are registered trademarks of 3D Systems, Inc.

**p/n 42-D136 Rev. A**