

Customer Information Bulletin



CIB# CIB00115

Date: 4/1/2019

Subject: New firmware update available – v1.0.4

NextDent 5100

There is a new firmware version available for the NextDent 5100 3D printer! It contains the latest updates from v1.0.2 and adds the **Resin-Tray Cleaning** function, which can be used after a failed build. Please watch the easy how-to video on this function at the link below:

[NextDent 5100](#)

To download/install the new firmware update:

Follow the procedure **Update Printer Firmware via USB** in the NextDent 5100 User Guide. The procedure can be found here: [NextDent 5100 Update Printer Firmware via USB](#).

DISCLAIMER OF LIABILITY: The following supersedes any related provision in your company's forms, letters, and agreements from, by or with 3D Systems Inc. 3D Systems Inc. makes no warranty, whether expressed or implied, including warranties of merchantability or of fitness for a particular purpose for this product. No statements or recommendations contained in the product literature are to be construed as inducements to infringe any relevant patent now or hereafter in existence. Under no circumstances shall 3D Systems Inc. be liable for incidental, consequential, special, or other damages from alleged negligence, breach of warranty, strict liability or any other theory, arising out of the manufacture, use, sale, or handling of this product. In no event shall the liability of 3D Systems Corporation for any claims arising out of the manufacture, use, handling, or sale of its products exceed an amount equal to the buyer's purchase price.

© Copyright 2019 by 3D Systems, Inc. All rights reserved. Subject to change without notice. The 3D logo is a registered trademark of 3D Systems.