

Customer Information Bulletin



CIB# CIB00132
Date: 6/18/2019

Subject: Introducing 3D Sprint 2.11

NextDent™ 5100

3D Sprint 2.11 for NextDent 5100 has now launched! Please review the document [What's New in 3D Sprint 2.11](#) to see all the latest features that this release includes. Some highlights for NextDent 5100 are below:

1. Build-time estimation has been improved.
2. Improved support structures have been added.
3. You can now check for and install firmware updates for your printers through 3D Sprint.
4. There are new and updated build styles for the following materials:
 1. Model 2.0 White Orthodontic – NEW!
 2. Tray – Shade Blue – UPDATED!
 3. Tray – Shade Pink – UPDATED!

DISCLAIMER OF LIABILITY: The following supersedes any related provision in your company's forms, letters, and agreements from, by or with 3D Systems Inc. 3D Systems Inc. makes no warranty, whether expressed or implied, including warranties of merchantability or of fitness for a particular purpose for this product. No statements or recommendations contained in the product literature are to be construed as inducements to infringe any relevant patent now or hereafter in existence. Under no circumstances shall 3D Systems Inc. be liable for incidental, consequential, special, or other damages from alleged negligence, breach of warranty, strict liability or any other theory, arising out of the manufacture, use, sale, or handling of this product. In no event shall the liability of 3D Systems Corporation for any claims arising out of the manufacture, use, handling, or sale of its products exceed an amount equal to the buyer's purchase price.

© Copyright 2018 by 3D Systems, Inc. All rights reserved. Subject to change without notice. NextDent 5100 is a trademark and the 3D logo is a registered trademark of 3D Systems.