

Customer Information Bulletin



CIB# CIB00192
Date: 5/11/2020

Subject: Introducing 3D Sprint 2.13

NextDent® 5100

3D Sprint 2.13 for Figure 4 Standalone and Figure 4 Jewelry has now launched! Please review the document [What's New in 3D Sprint 2.13](#) to see all the latest features that this release includes. The following are some of the key updates from this release:

- 1. Create Bars tool has been added** – This new tool allows you to quickly add an enforcement bar to horseshoe models and denture bases with a few simple clicks. Enforcement bars can help control material shrinkage and part deformation.
- 2. Stacked arches for Model 2.0 White** – This new tool is a one-click solution for automated preparation and placement of stacked parts. This optimized workflow enables printing 30 orthodontic dental models in under two hours with the new Ortho Stacked build style.
- 3. Auto update mechanism for new or updated materials** – Previously, new/updated materials and build styles had to be added to 3D Sprint via a separate update that was downloaded from Infocenter. Now, these updates will be added automatically upon availability.

3D Sprint 2.13 requires your printer to be running firmware v1.0.4-9. Please ensure you have updated your firmware to the latest version.

DISCLAIMER OF LIABILITY: The following supersedes any related provision in your company's forms, letters, and agreements from, by or with 3D Systems Inc. 3D Systems Inc. makes no warranty, whether expressed or implied, including warranties of merchantability or of fitness for a particular purpose for this product. No statements or recommendations contained in the product literature are to be construed as inducements to infringe any relevant patent now or hereafter in existence. Under no circumstances shall 3D Systems Inc. be liable for incidental, consequential, special, or other damages from alleged negligence, breach of warranty, strict liability or any other theory, arising out of the manufacture, use, sale, or handling of this product. In no event shall the liability of 3D Systems Corporation for any claims arising out of the manufacture, use, handling, or sale of its products exceed an amount equal to the buyer's purchase price.

© Copyright 2018 by 3D Systems, Inc. All rights reserved. Subject to change without notice. The 3D logo and Figure 4 are registered trademarks of 3D Systems.